


μFOCUS 450

THE MOST VERSATILE X-RAY SOURCE
FOR CUTTING EDGE XPS, HAXPES AND NAP-XPS

KEY FEATURES

- Multiple photon energies at 1487 eV and 2984 eV
- Chromium extension for HAXPES at 5414 eV
- Variable spot size from <math><100\ \mu\text{m}</math> to $>1000\ \mu\text{m}$
- Ultimate photon flux densities
- Full software control



SPECS[™]
A member of SPECSGROUP

μFOCUS 450

THE MOST VERSATILE X-RAY SOURCE
FOR CUTTING EDGE XPS, HAXPES AND NAP-XPS

Multi Wavelength Excitation

The μFOCUS 450 is a most versatile small-spot multi-wavelength X-ray monochromator for surface analysis and depth profiling applications. This fully computer controlled device allows for in situ switching between different emitters for Al, Ag and Cr excitation, thus providing X-rays with high flux density for various application such as small, medium and high energy XPS and surface layer thickness characterization.


Design

The new monochromator is an all in one solution for XPS and HAXPES systems, without the space consuming solutions of several individual light sources. It combines the performance of the μFOCUS 500/600, for Al with 1487 eV and Ag with 2984 eV excitation, and the μFOCUS 730, for Cr with 5414 eV excitation. The new design hosts up to three anodes and the corresponding monochromator optics in one Rowland Circle based housing. Switching between the anodes and the optics is fully computer controlled. The chromi-


um extension is hereby an additional add-on to the base version. With a variable spot size of less than 100 μm to more than 1 mm, this X-ray source fits most requirements, starting with small spot (and NAP) hemispherical analyzers to large or X-ray sensitive samples. A special NAP extension for operation at elevated pressures is available and ideally fits the new AEOLUS 150 NAP analyzer.

Laboratory XPS, HAXPES and NAP XPS

A reliable and functional X-ray source is the key to any photo electron spectrometer, as the best analyzer can't perform without a strong partner. The μFOCUS 450 is the solution for modern laboratories, either dedicated to complex and custom surface analysis, high energy excitation studies, modern material analysis under ambient conditions – or even all in one system. Its modularity allows the user to choose the correct configuration without compromising performance or future upgrades.



Two layer thin film on Si-wafer, depth profile gate oxide



Monochromator layout with the Rowland Circle; in mm


Depth Profiling

Non destructive analysis of layer thickness requires high energy excitation to overcome the inelastic mean free path and emit electrons from deeper layers for analysis. On the other hand, lower cross section and thus count rate reduce the efficiency of the overall experiment. Combining several excitation lines, from low to high energy, give a comprehensive and reliable result, which can efficiently be used for layer thickness analysis. Well focusing of the source enhances the flux density and hence, compensates the reduced efficiency.

Technical Data

Specification	Al K _α	Ag L _α	Cr K _α
Photon Energy	1487 eV	2984 eV	5414 eV
Minimum Spot Size	< 100 μm	< 100 μm	< 100 μm
Maximum Spot Size	> 1000 μm	> 1000 μm	> 1000 μm
Maximum Power	400 W	400 W	400 W (30 kV)
Max Flux (Ph/(s×mm ²))	1.9×10 ¹⁰	1.2×10 ⁹	1.9×10 ¹⁰
X-ray Line width (FWHM)	< 220 meV	< 450 meV	< 500 meV
Rowland Circle	450 mm	416 mm	676 mm

Dimensions



μFOCUS with chromium extension; dimensions in mm



NAP XPS

Due to signal absorption in the surrounding gas, NAP XPS is one of the most challenging applications in surface analysis. High performance analyzers, such as the AEOLOS 150 NAP, require a matched X-ray source, fitting the small analyzer field of view and the need for a high flux density on smallest spots. The μFOCUS 450 is designed to cover these requirements for highly efficient NAP XPS measurements.

SPECS Surface Nano Analysis GmbH
Voltastrasse 5
13355 Berlin / Germany
www.specs-group.com

T +49 30 46 78 24 0
F +49 30 46 42 083
E info@specs.com

Version 04.23

SPECSGROUP